



Great Salt Lake Water Quality Studies

Avian Blood Sample Analysis



Avian Blood Sample Analysis

Background

**Observed high Se concentrations in
nesting bird blood samples**

Verified methods and checked data

**Wanted to complete inter-lab comparison
of 2006 method to verify results**



Avian Blood Sample Analysis

**No surplus blood samples from
summer 2006 nesting season**

**Therefore no samples available to
split with USGS to verify Se blood
levels**



Avian Blood Sample Analysis

56 Eared Grebes are awaiting analysis

- 25 from early collection (planned for 20)
- 31 from late collection (planned for 20)

**Adequate blood volume from each bird
for Se analysis but not enough for
inter-lab comparison and additional
analyses**

**Therefore need to composite some blood
samples**



Avian Blood Sample Analysis

**Mike Conover and Gary Santolo
selected “spare” samples to be
used for compositing**

**Tom May/USGS will create 3 freeze-
dried composites with 5 blood
samples in each, split samples to
be sent back to LET**



Avian Blood Sample Analysis

USGS and LET to complete the following inter-lab comparisons:

- Current HGAA method used for Se blood analysis
- Proposed method for Se + Hg analysis
- Standard Reference Materials for Se/Hg



Avian Blood Sample Analysis

Se + Hg Method for Analysis

- 10 eared grebes from early season and 10 from late season from the Hat Island sampling location
- 10 goldeneyes from late season
- Undetermined number of nesting birds from May/June 2007



Avian Blood Sample Analysis

Se + Hg Method for Analysis

For LET to be able to complete Se + Hg on eared grebe blood samples, LET must change prep procedure

Current Prep Procedure:

- Freeze-dry samples
- Perform dry-ash digestion for Se analysis



Avian Blood Sample Analysis

Se + Hg Method for Analysis

Proposed Prep Procedure:

- Freeze-dry samples
- Perform microwave digestion using nitric/peroxide solution
- Split digestate for Hg and Se analysis
- Perform dry-ash digestion for Se analysis



Avian Blood Sample Analysis

LET and USGS will compare results for existing Se method and new Se + Hg prep method

After inter-lab comparison is complete and satisfactory, LET will:

- 1. Analyze 20 grebe blood samples for Se (with existing prep method)**
- 2. Analyze 20 grebe blood samples for Se + Hg using new prep method**



Avian Blood Sample Analysis

Next Steps:

1. Identify nesting species and number of birds to be collected/analyzed for Se + Hg
2. Determine path forward



Avian Blood Sample Analysis

Next Steps:

1. If we confirm high Se in blood, then what?
2. If we have low Se in blood, then what?
3. If we confirm high Se and high Hg in blood, then what?
4. If we confirm high Se but low Hg in blood, then what?